



XRF SAMPLE CUP “PALM-HELD” SNAP-ON RING and SLEEVE FASTENER

HOME
BACK TO SNAP-ON RING FASTENER
VENT HOLE PUNCH
TRANSFER PIPETTES

Held in the palm, this handy device effortlessly aids in attaching “Snap-On Rings” and “Sleeves” to XRF Sample Cups. Averts contamination by contact during assembly, allows uniform application of pressure, eliminates Snap-On Ring and Sleeve misalignment and ensures taut thin-film sample support planes; ideally suitable for multiple sample preparations.



Usable With the Following Listed XRF Sample Cups

Sample Cup Cat. No.	Sample Cup Cat. No.	Sample Cup Cat. No.	Sample Cup Cat. No.
1330	1530	1830 SE	1060
1330 SE	1530 SE	1840	1065
1340	1540	1850	1070
1430	1730	1930	1075
1430 SE	1730 SE	1930 SE	1080
1440	1740	1940	1083
1440L	1830	1940L	1085

SPECIFICATIONS AND ORDERING INFORMATION

CATALOG NUMBER	DESCRIPTION
1325	Palm-Held Snap-On Ring and Sleeve Fastener

Chemplex[®] INDUSTRIES, INC.

2820 SW 42nd Avenue
Palm City, Fl. 34990-5573, USA

Continental USA Toll Free: (800) 4-CHEMPLEX
Phone: (772) 283-2700
Fax: (772) 283-2774

Email: sales@chemplex.com

On Line Sales: <http://www.chemplex.com>