



XRF SAMPLE CUP “SNAP-ON RING” FASTENER

HOME
PALM HELD FASTENER
VENT HOLE PUNCH
TRANSFER PIPETTES

This unique device is designed to facilitate the attachment of thin-film sample support windows to sample cups utilizing “Snap-On Rings.” Averts contamination by contact during assembly, allows uniform application of pressure, eliminates Snap-On Ring misalignment and ensures taut thin-film sample support planes; ideally suitable for multiple sample preparations.



SIMPLE, FAST AND CONVENIENT TO USE



1. Insert sample cup on platform
2. Cover sample cup opening with thin-film
3. Position “Snap-On Ring” over thin-film
4. Depress lever and release



USABLE WITH THE FOLLOWING XRF SAMPLE CUPS

SAMPLE CUP CAT. NO.	SAMPLE CUP CAT. NO.
1330	1730 SE
1330 SE	1740
1430	1830
1430 SE	1830 SE
1440	1840
1530	1930
1530 SE	1930 SE
1540	1940
1730	1850

SPECIFICATIONS AND ORDERING INFORMATION

CATALOG NUMBER	DESCRIPTION
1300	Snap-On Ring Fastener

Chemplex[®] INDUSTRIES, INC.

2820 SW 42nd Avenue
Palm City, Fl. 34990-5573, USA

Continental USA Toll Free: (800) 4-CHEMPLEX
Phone: (772) 283-2700
Fax: (772) 283-2774

Email: sales@chemplex.com

On Line Sales: <http://www.chemplex.com>